



**Wellesley  
Cosmetic  
Surgery**

Arthur Shektman, MD  
332 Washington St.  
Suite 340  
Wellesley MA 02481  
Tel: 781-239-0680  
Fax: 781-239-0681

Email: doctor@  
wellesleycosmeticsurgery.com

## Drain Instructions

After your procedure you will find a plastic tube extending from under your garment and attaching to a plastic bulb. Initially the fluid coming out will be red, but after a day or two it will turn yellow. When the bulb fills up you will need to empty it. It is important to do it properly, or the drain will not work!

On the top of the bulb you will find a small nozzle closed off with a flap.

1. Open the flap
2. Turn the bulb over and squeeze the fluid out of it
3. Continue to squeeze the bulb so that all air is expelled from it
4. While keeping the bulb squeezed replace the flap and close off the nozzle
5. Make sure the flap is sealed tightly over the nozzle

If you did it correctly the bulb will remain collapsed, which means that there is vacuum inside and the suction will work properly



Certified by the American  
Board of Plastic Surgery

Facility accredited by the Joint  
Commission on Accreditation  
of Healthcare Organizations